

## LIFE SCIENCES

Nanoscribe's Photonic Professional GT 3D printers allow to easily produce tailored and reproducible 3D scaffolds and matrices. Here, the tailoring is achieved by on one hand a wide range of processable materials and on the other hand a resolution and feature sizes which, in combination, allow to mimic in vivo physiological environments. Accessible patterning scales cover sub-micrometer features and sample volumes of up to several 10 mm<sup>3</sup> at the same time. Samples can be produced from biocompatible, cell binding/repelling, hydrophobic and hydrophilic polymers with different mechanical characteristics as well as hydrogels.




### CELL / EXTRACELLULAR MATRIX INTERACTION

**Challenge:** Fabrication of a microstructured 3D cell scaffold consisting of a cell-repellent material (grey) and a cell adhesion promoting material (red) as per design.

**Solution:** First, a core scaffold made of PEG-DA mixed with PETA (grey) is additively printed. Then, a second printing step adds cell adhesion points made of OrmoComp (red) precisely at the designed locations.

**Source:** DOI: 10.1002/adma.201004060




### TARGETED DRUG DELIVERY

**Challenge:** One-step fabrication of cargo-loaded helical swimming microrobots made of biodegradable superparamagnetic hydrogel for minimal invasive intervention.

**Solution:** A hydrogel consisting of PEG-DA/PETA, magnetite nanoparticles and photoinitiator is 3D printed. For loading, the helices are immersed in methylene blue.

**Source:** DOI: 10.1002/adma.201503112




### CELL MIGRATION

**Challenge:** Fabrication of cage-like 3D micro-structures with sufficiently small pore size to differentiate tumorigenic from non-tumorigenic human breast cancer cells based on different invasion behavior.

**Solution:** Cubic cages with meshes of approximately 3, 18 and 85 µm<sup>2</sup> pore size as side faces are printed in Nanoscribe's biocompatible photoresist IP-L 780.

**Source:** DOI: 10.1038/srep10531




### TISSUE ENGINEERING / 3D CELL CULTURE

**Challenge:** High-fidelity replication of a human trabecular bone structure from high resolution 3D data obtained from a µ-CT scan.

**Solution:** µ-CT data is transferred to the standard 3D printing workflow. The scanned 3D object is reproduced with sub-micrometer precision in OrmoComp.

**Source:** DOI: 10.1016/j.actbio.2014.05.032



Learn more about the world's  
highest resolution 3D printer  
Photonic Professional 

### Nanoscribe GmbH


Phone +49 721 981 980 0

E-Mail lifesciences@nanoscribe.com

Web www.nanoscribe.com

## MECHANICAL METAMATERIALS

Mechanical properties of materials can be tailored by material composition but also by geometrical design. The field of mechanical metamaterials aims at exploiting exactly this twofold design freedom. Nanoscribe's Photonic Professional GT offers three-dimensionality (3D), submicron feature size and resolution as well as hierarchical structuring capability. Auxetics and mechanical cloaks are successfully demonstrated by polymeric structures. By subsequent post-processes such as surface coating or electro-plating, damage-tolerant lightweight materials are also made of composites, metals and ceramics.




### CHIRAL MECHANICAL METAMATERIAL

**Challenge:** Examine non-Cauchy elastic metamaterials that can respond to compression with a twisting motion (instead of buckling/bending).

**Solution:** Scaling behavior is studied on several, easily printed samples with identical overall dimensions but consisting of different numbers of unit cells.

**Source:** DOI: 10.1126/science.aao4640




### HIGH-STRENGTH CERAMIC POLYMER COMPOSITE

**Challenge:** Fabrication of a 3D cellular material with low filling factor built of submicron elements demonstrating size-dependent strengthening effects.

**Solution:** 3D polymeric structures with submicron feature size are printed. Ceramic composites are then obtained by subsequent ALD coating with alumina.

**Source:** DOI: 10.1073/pnas.1315147111




### RESILIENT 3D HIERARCHICAL METAMATERIAL

**Challenge:** Manufacturing of 3D fractal-like architectures with features in the sub-micron regime for the attainment of various mechanical properties, e.g. ultralightweight, recoverability.

**Solution:** 3D hierarchical nanolattices with individual beams comprised of multiple self-similar unit cells with length scales spanning over four orders of magnitude are fabricated.

**Source:** DOI: 10.1073/pnas.1509120112




### AUXETICS / NEGATIVE POISSON RATIO MATERIAL

**Challenge:** Micro-manufacturing of a material that laterally expands when stretched and shrinks when compressed, for applications in biomimetics or magneto-mechanical microsystems.

**Solution:** 3D complex geometries with submicron details are fabricated precisely, reproducibly and mechanically stably.

**Source:** DOI: 10.1088/0964-1726/23/8/085033



Learn more about the world's  
highest resolution 3D printer  
Photonic Professional 

#### Nanoscribe GmbH

Phone +49 721 981 980 0

E-Mail [m.metamaterials@nanoscribe.com](mailto:m.metamaterials@nanoscribe.com)

Web [www.nanoscribe.com](http://www.nanoscribe.com)

## MICROFLUIDICS

The precise control and manipulation of liquids in very small volumes is the essence of microfluidics, for example in lab-on-a-chip applications. Nanoscribe's Photonic Professional *GT* offers the highest resolution available in 3D printing today, enabling tailored and smooth surfaces to be producible as desired. The additive fabrication approach overcomes technical challenges such as high aspect ratios or high surface to volume ratios. Thus, 3D intertwined compact mixers and filter elements can be achieved. By fully exploiting the design freedom even biomimetic micro-needles for painless drug delivery and customized nozzles can be printed.




### 3D MICROFLUIDIC FILTER

**Challenge:** Fabrication of a non-clogging micro-porous 3D filter with connectors in thick-film SU-8.

**Solution:** The Photonic Professional *GT* allows for high resolution 3D printing over the total height of the filter element.

**Source:** [www.nanoscribe.de/en/applications/microfluidics/](http://www.nanoscribe.de/en/applications/microfluidics/)




### TAILORED SURFACES WITHIN MICROFLUIDIC DEVICES

**Challenge:** Introduction of micro-textured surfaces into microfluidic channels to study impact on cell motility or to design hydrophobic/hydrophilic characteristics.

**Solution:** 3D printing allows for tailor-made sub-micrometer modifications of the surface topography. For further replication, these can be used as masters.

**Source:** DOI: 10.1007/s10544-014-9864-2




### MASTER OF FILTER ELEMENT

**Challenge:** Rapid and precise fabrication of masters for injection molding or imprinting with heights and aspects ratios exceeding those accessible with other techniques.

**Solution:** Additive fabrication easily overcomes the design constraints while enabling rapid design iterations and preserving the replicability for mass production.

**Source:** [www.nanoscribe.de/en/applications/microfluidics/](http://www.nanoscribe.de/en/applications/microfluidics/)




### NOZZLES

**Challenge:** Fabrication of designer microfluidic nozzles for the precise formation of droplets, sprays or jets.

**Solution:** Almost arbitrarily shaped nozzles with micrometer precision can be printed. This allows for internal features and precisely structured orifices.

**Source:** DOI: 10.1364/OE.24.011515



Learn more about the world's  
highest resolution 3D printer  
Photonic Professional *GT*

### Nanoscribe GmbH

Phone +49 721 981 980 0

E-Mail [microfluidics@nanoscribe.com](mailto:microfluidics@nanoscribe.com)

Web [www.nanoscribe.com](http://www.nanoscribe.com)

## MICRO-OPTICS

Nanoscribe's 3D printer Photonic Professional *GT* allows producing almost arbitrary micro-shapes with optically smooth surfaces in an additive and tool-free way. This effectively circumvents limitations imposed by mechanical tools, and geometrical or process design-constraints often encountered with techniques such as subtractive machining, (greyscale) lithography, photoresist reflow and wet-etching. In this way, steep slopes for high numerical aperture micro-lenses, arrays with high filling factors and varying curvatures, as well as more complex 3D shapes can be achieved. Mass replication is possible by using metal replicas of these shapes.




### VERTICAL SIDEWALLS

**Challenge:** Hemispherical micro lenses with steep slopes to achieve low f-numbers in a closely spaced array.

**Solution:** The layer by layer writing process allows for smooth surfaces over a 0° to 90° angular range, and even for negative slopes (undercuts).

**Source:** [www.nanoscribe.de/en/applications/micro-optics/](http://www.nanoscribe.de/en/applications/micro-optics/)




### HIGH FILLING FACTOR / COMPACT ARRAY

**Challenge:** Microlens arrays with increased optical efficiency due to closely packed lenses.

**Solution:** 3D printing allows overlapping lenses, enabling designs with 100% filling factor. Arbitrary grids are possible and the lens shape can be varied throughout the array.

**Source:** [www.nanoscribe.de/en/applications/micro-optics/](http://www.nanoscribe.de/en/applications/micro-optics/)




### SHARP EDGES AND INTERNAL ANGLES

**Challenge:** Fabricate corner cube arrays on the microscopic scale with well-defined edges.

**Solution:** The absence of a mechanical tool and the highly localized polymerization from two-photon absorption allow for internal angles and the small feature sizes necessary for sharp corners.

**Source:** [www.nanoscribe.de/en/applications/micro-optics/](http://www.nanoscribe.de/en/applications/micro-optics/)




### MASS REPLICATION VIA NI-SHIM

**Challenge:** Fast and low cost production of micro-optical elements.

**Solution:** A nickel shim can be fabricated from the printed polymer structures by electroforming, allowing standard replication techniques such as injection molding or hot embossing to be used.

**Source:** [www.nanoscribe.de/en/applications/micro-optics/](http://www.nanoscribe.de/en/applications/micro-optics/)



Learn more about the world's  
highest resolution 3D printer  
Photonic Professional *GT*

#### Nanoscribe GmbH


Phone +49 721 981 980 0

E-Mail [microoptics@nanoscribe.com](mailto:microoptics@nanoscribe.com)

Web [www.nanoscribe.com](http://www.nanoscribe.com)

## MICROBOTICS

Nanoscribe's 3D microprinter Photonic Professional *GT* allows for microfabrication with an intrinsic capacity to build complex and custom 3D micro tools. It enables the fabrication of microrobots that interact with tissue and single living cells and conduct medical tasks at the micro scale. Remotely controlled micromachines, which can be powered by magnetic, chemical or optical actuation, are crucial for minimal invasive operations. Depending on the desired actuation mechanism microrobots can be printed from photopolymers, nanoparticle composites as well as hydrogels. Metallic coatings can be applied if required.




### MICROTRANSPORTERS FOR TARGETED DRUG DELIVERY

**Challenge:** Controlled loading, transport and release of therapeutic agents in fluids by means of a microtransporter remotely actuated by a magnetic field that rotates a magnetic shaft.

**Solution:** The complete transporter is printed without further assembly. A sacrificial printed shell allows for selective coating of the structure with Ni/Ti.

**Source:** DOI: 10.1002/adma.201503095




### FERTILIZATION-ASSISTING MICROMOTORS

**Challenge:** Cell-compatible micromachine able to catch, transport and release living sperm to an oocyte overcoming low sperm motility-induced fertility issues.

**Solution:** Printed polymer microhelices are coated with nickel and titanium allowing the helices to act as motors in a magnetic field and providing biocompatibility, respectively.

**Source:** DOI: 10.1021/acs.nanolett.5b04221




### OPTICAL MICROSYPHINGE

**Challenge:** Micro tool capable of movement with six degrees of freedom and controlled as well as precise cargo loading, transport and release.

**Solution:** A hollow micro vessel with several trapping handles is printed and moved by optical tweezing. Optically changing the pressure inside the vessel allows for loading and unloading of particles.

**Source:** DOI: 10.1038/lisa.2016.148




### LIQUID CRYSTAL ELASTOMER (LCE) MICROWALKER

**Challenge:** Soft microrobot that reversibly deforms and is powered by light for walking and jumping, overcoming strong surface adhesion.

**Solution:** Conical legs are printed on a previously printed LCE structure to reduce adhesion. Modulated laser light contracts the robot body, changing the tilt of the legs leading to locomotion.

**Source:** DOI: 10.1002/adma.201501446



Learn more about the world's  
highest resolution 3D printer  
Photonic Professional *GT*

### Nanoscribe GmbH

Phone +49 721 981 980 0

E-Mail [microrobotics@nanoscribe.com](mailto:microrobotics@nanoscribe.com)

Web [www.nanoscribe.com](http://www.nanoscribe.com)

## 3D PHOTONIC MATERIALS

Nowadays, in telecommunication technology as well as for photonic circuits and devices, materials whose optical response can be tailored by design are being developed. Nanoscribe's Photonic Professional *GT* offers the highest resolution known in 3D fabrication and submicron feature sizes typically on or below the length scale of the wavelength of light, leading to novel optical properties proven in reflection and transmission measurements. Subsequent processes such as electroplating, and single or double inversion using ALD and CVD allow to transfer polymer templates into different materials, e.g. gold, titanium dioxide, silicon.




### 3D INVISIBILITY CLOAK

**Challenge:** Realization of spatially varying, tailored effective refractive index distribution in 3D for hiding a metallic bump at optical wavelengths.

**Solution:** 3D woodpile photonic crystal with rod distance of 800 nm and varying rod width is fabricated with the polymer material IP-L.

**Source:** DOI: 10.1126/science.1186351




### SILICON PHOTONIC BANDGAP MATERIAL

**Challenge:** Rationally designed aperiodic 3D structure from a high refractive index material with a photonic bandgap in the near infrared.

**Solution:** Polymer networks are printed, coated with TiO<sub>2</sub> by ALD and the polymer is removed by calcination. The remaining shrunk TiO<sub>2</sub> structure is infiltrated with silicon by CVD.

**Source:** DOI: 10.1002/adom.201300415




### PHOTONIC PARTICLE-ACCELERATOR

**Challenge:** Fabrication of an air waveguide inside a 3D photonic-crystal in order to accelerate electrons using a confined optical mode.

**Solution:** Woodpile structure with a waveguide inside is printed. After fabrication, the silicon replica of the printed structure is generated via silicon double inversion.

**Source:** DOI: 10.1364/OE.20.005607




### TAPERED GOLD-HELIX BROADBAND CIRCULAR POLARIZER

**Challenge:** Fabrication of arrays of gold-helices whose radii increase gradually along the helix axis serving as broadband circular polarizers in the mid-infrared.

**Solution:** Positive-tone resist template with helices is printed and then transferred to gold by LIGA processes. The resist is then removed by plasma etching.

**Source:** DOI: 10.1063/1.3693181



Learn more about the world's  
highest resolution 3D printer  
Photonic Professional *GT*

#### Nanoscribe GmbH

Phone +49 721 981 980 0

E-Mail [photonicmaterials@nanoscribe.com](mailto:photonicmaterials@nanoscribe.com)

Web [www.nanoscribe.com](http://www.nanoscribe.com)

## PLASMONICS

Nanoscribe's Photonic Professional *GT* system is the world's highest resolution 3D printer. It offers resolution between what can be achieved with electron beam lithography and UV lithography, providing at the same time the full fabrication flexibility for 2D, 2.5D as well as for 3D objects. Highest accuracy for high performance demanding researchers is enabled with one single instrument that can be installed in a non-clean room as well as in a clean room environment. Maskless lithography is possible on silicon wafers as well as on non-conductive substrates like glass. The four examples below show ways to achieve plasmonic behavior.




### PERFECT ABSORBER

**Challenge:** Fabrication of a metamaterial absorber in the MID-IR range, that shows good heat dissipation.

**Solution:** Printing of a polymeric four-tined fish-spear-like resonator array that is covered by a continuous metallic film for excellent thermal conduction.

**Source:** DOI: 10.1002/adma.201300223



Courtesy of Prof. Mu Wang, Nanjing University, China

### PLASMONIC COLOR DISPLAY

**Challenge:** Develop a reflective display with pixels based on a surface plasmon resonance that can be tuned electrically to modulate the color.

**Solution:** Fabrication of a 2.5D master pattern array of nanowells and subsequent processes like nanoimprinting, aluminum coating and placing of liquid crystals.

**Source:** DOI: 10.1038/ncomms8337





Image: D. Franklin et al., Nature Communications 6, 7337 (2015)

### SUBWAVELENGTH LIGHT FOCUSING

**Challenge:** Provide a novel methodology for realization of 3D plasmonic focusing with a signal stronger than other techniques for subwavelength light focusing.

**Solution:** Fabrication of hollow pyramids with apertures of varying sizes and subsequent coating with a thin gold film on different substrates including a copper grid.

**Source:** DOI: 10.1364/OE.23.022564




Courtesy of Prof. Changzhi Gu, IOP, CAS

### PLASMONIC NANOANTENNA

**Challenge:** Establish an easy-to-use alternative for electron beam lithography for fabrication of plasmonic nanoantennas with line widths down to 105 nm.

**Solution:** The Photonic Professional *GT* enables patterning of polymer serving as an etch mask on a gold layer covering an area of tens of mm<sup>2</sup>.

**Source:** DOI: 10.1021/acsp Photonics.5b00141



Courtesy of 4th Physics Institute, University of Stuttgart



Learn more about the world's  
highest resolution 3D printer  
Photonic Professional *GT*

### Nanoscribe GmbH


Phone +49 721 981 980 0

E-Mail [plasmonics@nanoscribe.com](mailto:plasmonics@nanoscribe.com)

Web [www.nanoscribe.com](http://www.nanoscribe.com)

## BIOMEDICAL ENGINEERING

Nanoscribe's Photonic Professional *GT* 3D printers offer full design freedom and submicron precision for the fabrication of customized biomedical devices. The high resolution of the printers is especially favorable for mimicking natural environments at sub-cellular scale for lab-on-a-chip applications. Moreover, the printer is used to develop precision instruments needed in minimally invasive operations. Biocompatibility of multiple 3D-printed parts has been demonstrated in contact with living cells and tissue. The printable materials range from photoresins to biodegradable hydrogels and functional resins containing nanoparticles.




### BLOOD-BRAIN BARRIER MODEL

**Challenge:** Fabricating a real-scale biohybrid model of the blood-brain barrier for drug screening. This model requires tubular structures of 10  $\mu\text{m}$  diameter and pores of 1  $\mu\text{m}$  diameter as supporting scaffold for endothelial cells.

**Solution:** A millimeter-long system of porous cylindrical channels is 3D printed. The cells cultivated on this scaffold build a biological barrier mimicking the blood-brain barrier.

**Source:** DOI: 10.1002/sml.201702959




### NERVE INTERFACE IMPLANTS

**Challenge:** Producing micro-scale implants to stimulate and track nerve responses.

**Solution:** Nanoclips with trapdoors for fixing nerves, slots for placing electrodes and through-holes for surgical handling are 3D printed. Transplanted into zebra finches they show healthy nerve activity post implant.

**Source:** DOI: 10.1088/1741-2552/aa5a5b




### MULTI-LENS OPTICS ON A FIBER

**Challenge:** Fabrication of a compound lens system consisting of three objectives stacked on the tip of 125  $\mu\text{m}$  diameter optical fiber. The device can find applications in endoscopy.

**Solution:** A triplet lens system is fabricated in one printing step without the need of further assembly nor alignment of the lenses. Printing is done in situ on top of the fiber.

**Source:** DOI: 10.1038/nphoton.2016.121




### REPLICATED HOLLOW MICRONEEDLES

**Challenge:** Replication of hollow microneedle arrays with microfluidic channels and sharp tips. The replicas aim to be used for point-of-care collection of blood or drug delivery.

**Solution:** Needles are additively manufactured from a CAD model. The needles' master is cast in a negative mold. The mold is replicated by soft embossing into a thermoplastic cycloolefin polymer that shows agent delivery into a rabbit's ear.

**Source:** DOI: 10.1038/micronano.2017.34



Learn more about the world's  
highest resolution 3D printer  
Photonic Professional *GT*

#### Nanoscribe GmbH

Phone +49 721 981 980 0


E-Mail [biomedical@nanoscribe.com](mailto:biomedical@nanoscribe.com)

Web [www.nanoscribe.com](http://www.nanoscribe.com)



## BIOMIMETICS

Natural systems with millions of years of evolution inspire technological innovation in the design and development of functional and smart materials. Nanoscribe's Photonic Professional *GT* offers highest resolution in 3D printing and geometric design freedom to emulate biological structures at nano-, micro- and millimeter scales. Shape, size and distribution of the structures can be easily changed in order to investigate and optimize materials. Due to structural design at the microscale, materials can adopt properties such as hydrophobicity, elasticity or coloration – mimicking intricate architectures and finest features observed in nature.




### HIGHLY HYDROPHOBIC SURFACES

**Challenge:** 3D micropattern capable of air retention, hydrophobicity and water condensation required for fog collection inspired by salvinia molesta leaves.

**Solution:** Periodic arrays of crown-like hairs, hundred times smaller than the natural model, are printed. Structures are stable when submerged in water and air trapping is demonstrated.

**Source:** DOI: 10.1021/acsami.5b07722




### NONIRIDESCENT PHOTONIC NANOSTRUCTURES

**Challenge:** Tarantula hair-like nanostructured material which exhibits angle independent reflectance spectra due to microscale rotational symmetry and hierarchy.

**Solution:** Hierarchical multilayer cylinders are fabricated out of a transparent polymer material with layer thicknesses of 300 nm and layer spacing of 450 nm that show a green hue up to 70° incidence angle.

**Source:** DOI: 10.1002/adom.201600599




### MICRONEEDLES FOR PAINLESS BLOOD COLLECTION

**Challenge:** Microneedles mimicking the mosquito's proboscis for painless and easy insertion into the skin without bending or collapsing due to reduced puncturing resistance.

**Solution:** Microneedle halves of mm-length are printed with sharp tips, spikes and small holes. Two assembled halves advance alternately and can penetrate artificial skin to suck blood.

**Source:** DOI: 10.1109/TRANSDUCERS.2015.7180876




### CELL-ACTUATED HYDROGEL MICROSKELETON

**Challenge:** Stable and deformable hydrogel scaffold that allows selective cell adhesion and differentiation of C2C12 cells into mature muscle fibers that contract and relax when voltage is applied.

**Solution:** The printed microskeleton is made of a PEGDA-based hydrogel mixed with PETA which promotes stiffness and protein adhesion. The cell-seeded bio-MEMS system deforms reversibly under electrical stimulus.

**Source:** DOI: 10.1109/MEMSYS.2017.7863350



Learn more about the world's  
highest resolution 3D printer  
Photonic Professional *GT*

### Nanoscribe GmbH

Phone +49 721 981 980 0

E-Mail biomimetics@nanoscribe.com

Web www.nanoscribe.com